

# Original Manufacturer Chemistry Reagents



## Modular System

Hepatic	Renal	Cardiac	Inorganic &Anemia	Lipid	Immune protein	Diabetes	Rheumatism	Pancreatitis	Lung
ALT	CREA	CK	Mg	TC	C3	Glu	RF	AMY	ADA
AST	UREA	CK-MB	P	TG	C4	HbA1C	ASO	LIP	ACE
GGT	UA	a-HBDH	Ca	HDL-C	IgA	β-HB	CRP		
ALP	CO2	LDH	FER	LDL-C	IgG	FUN			
ALB	Cys-C	HS-CRP	Fe	Apo-A1	IgM				
TP	β2-MG	MYO	UIBC	Apo-B	IgE				
PA	RBP	D-Dimer	TRF	Lp(a)					
TBA	MALB		G6PD						
CHE	TPUC								
BIL-T									
BIL-D									
AFU									

High Performance & Unique Tests

HbA1C


ADA (American Diabetes Association) Endorses HbA1c for Diabetes Diagnosis




Advantages

- Enzymatic method
- Wide measuring range:3~16%
- Excellent sensitivity: Hb (15 μmol/L) , HbA1C (3 umol/L)
- Excellent precision with intra CV value < 2.0%
- Traceability to IFCC/NGSP/JCCL reference methods

Precision Contrast





The bias value between the enzymatic HbA1c test results and the stated value : -0.2 ~ 0.4%




The bias value between the immunoassay HbA1c test results and the stated value : -0.5 ~ 0.6%

Method comparison



Reagent contamination



- Sampl①: 0.9% saline    Sampl②: hemolysis solution    Sampl③: hemolysis solution+ enzymatic HbA1c reagents  
Sampl④: hemolysis solution+ immunoassay reagent I    Sampl⑤: hemolysis solution+ immunoassay reagent II

International standardization

International standardization certificates of HbA1c from NGSP.  
More information refers to webstie.(<http://www.cdc.gov>)  
NGSP(National Glycohemoglobin Standardization Program)



# β-HB (β-Hydroxybutyrate )

ADA (American Diabetes Association)  
Endorses β-HB for diagnosing and monitoring ketoacidosis

## Enzymatic Colorimetric Method

Liquid ready-to-use reagent  
Enzymatic Colorimetric Method for the measurement of β-Hydroxybutyrate, making automation easier  
Wide measuring range with 0.03~5.5 mmol/L  
Excellent sensitivity with minimal detection level as 0.03mmol/L  
To test β-Hydroxybutyrate – the major ketone in the blood, making it a very reliable test  
Measurement of ketones in serum rather than in urine helps eliminate the risk of false negatives due to insensitivity and false positives due to drug interference

## Traditional diagnostic method

**Sodium nitroprusside reaction(urine strip or tablet)**  
Qualitative result  
False-positive results(Sulfhydryl drugs, including the ACE inhibitor captopril),false-negative results( High dosage of Vitamin C)  
Not early period diabetic ketoacidosis diagnosis, for not testing β-HB  
Misleads doctors for the increasing acetoacetate during treatment ( β-HB oxidized into acetoacetate)

# ADA ( Adenosine Deaminase)

An Ideal Biomarker For Tuberculosis Screening

## Advantages


High sensitivity: Measured Low limit of detection 1U/L  
Wide linearity range: Measured 1~200 U/L  
Excellent precision with low CV value.  
Strong anti-interfering ability for lipemia, bilirubin, hemaglobin and ascorbic acid  
Good correlation to reference method  
Calibrators standardized to International Reference Material ERM-AD455/IFCC



## Comparison ADA assay with prevalent method

### Interferon-γ Release Assays

Manual Operation  
16-24 hrs( but longer if run in batches) to get results




### Polymerase Chain Reaction(PCR)

Very expensive  
Complicated operation

### Adenosine deaminase Assay

Automatic chemistry method  
Rapid test  
Cost effective  
Simple operation



# G6PD ( Glucose-6-phosphate Dehydrogenase)


An important parameter for G6PD deficiency screening

## Advantages

- Liquid reagent, ready to use
- Wide linearity range: 0~3000 U/L.
- Excellent precision.
- Excellent on-board stability: 21 days.
- Long shelf life: 15 months.
- The sample type is the backlog red blood cells.
- There is no interfering for lipemia, bilirubin, drugs and so on.



## The comparison with prevalent methods



# AFU ( α-L-fucosidase )

A sensitive marker for diagnosing hepatocellular carcinoma (HCC)

## Advantages

- Liquid reagent, ready to use
- High sensitivity and wide linearity range
- Excellent precision with CV< 2%
- Excellent stability
- Strong anti-interfering ability of ascorbic acid and intra-lipid.
- Calibrators standardized to International Reference Material ERM-AD455/IFCC.



The combined detection of AFU and alpha-fetoprotein (AFP) is better than AFP individual detection

- The combination test of AFP & AFU enzyme increases the diagnostic accuracy to 90% inspite of AFP alone.
- About58% patients with HCC negative for AFP could be diagnosed.
- In the HCC patients with tumor size < 3 cm, the concentration of AFP is very low and some of the results were negative. AFU is a good marker than AFP in the detection of early HCC.

	Sensitivity	Accuracy
AFP	52%	74%
AFU	70%	78%
APF+AFU	80%	90%